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Model : Waters

I. Purpose:

1. Water-soluble polymer molecular weight distribution determination.
2. Separation of water-soluble polymer and calculate the absolute molecular weight.

II. Sample preparation instructions and precautions:

1. Mobile phase should be prepared with ultra-pure water.
2. The aqueous mobile phase and samples solution should be microfiltered with a 0.2 μm membrane filter and degassed.
3. Flow rate of mobile phase should not be higher than 0.8 ml / min.
4. No air bubble is allowed in the injected sample.
5. Injection needles, solvent and membrane filter should be provided by the user.
6. Suggested sample concentration: 0.5 ~ 2 mg / ml, injection volume 200 μl .

III. Restrictions:

1. Only for suitable water-soluble polymer samples.
2. Sample should be dissolved in ultra-pure water.

IV. Application:

1. Registration for application in advance is required. Users are required to operate the instrument themselves.
2. For users outside the department application should be signed by the chairman of the Chemical Engineering Department and the professor-in-charge for the equipment.

V. Professor-In-Charge: Yu-Cheng Chiu