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- (3) Large sample chamber with the optional ability to cover the incident and scattered beam paths with windows for samples in ambient or other non-vacuum environment. (大型的樣品室，具含蓋入射及散射光束路徑的選擇能力，並開有視窗，以利常壓或非真空環境下的樣品測量。)
- (4) The sample changer drive mechanism is mounted outside the sample chamber to maximize the available space inside the chamber. (樣品更換器驅動機制是安裝於樣品室之外，以使樣品室內的可用空間最大化。)
- (5) An x-ray photodiode is provided inside the sample chamber to be used for performance diagnostics and also fine tuning the alignment of the monochromator optic and the pinhole positions. (提供一個 X 射線光電二極管於樣品室內，用以診斷操作性能，及微調單光器光學系統是否對齊及針孔位置。)
- (6) Video display of the beam position on the sample simplifies the alignment of small samples. (打在樣品上的光束位置，可於影像上顯示，這簡化了小樣品的對齊調整過程。)
- (7) The beamstop is embedded with an x-ray sensitive photodiode for in-situ transmission measurement. (光束停止點嵌入了一個對 X 射線敏感之光電二極管，以提供穿透 X 光強度的即時測量功能。)

#### **Detector(偵測器):**

2D Multiwire, gas-filled proportional type with (二維、多接線、充氣比例類型具有):

- (1) 120mm active diameter (120 毫米之有效直徑。)
- (2) Sub-microsecond time stamping for performing time resolved or periodic pump/probe experiments without waiting for data storage. (次微秒時間印記功能，以執行時間分離或週期性之泵送／探測實驗，不需等待資料儲存。)

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- 室溫測量條件下，固體樣品試片或薄膜尺寸 (size suggested for solid specimen or film under room temperature)：3 to 5 cm (length), 0.5 to 1 cm (width), and less than 1 to 2 mm (thickness). (The sample will be fixed on a sample rack with 42 holes (about 1 cm in diameter) on it using a 3M Scotch tape.)
- 變溫測量條件下(20-300C)，固體樣品試片或薄膜尺寸 (size suggested for solid specimen or film under varied temperature (20-300 C))：0.5 cm (length), 0.5 cm (width), and less than 1 to 2 mm (thickness). (The sample will be placed on a temperature-controlled liquid sample cell with sample size limitation mentioned above.)
- 固體粉末數量 (amount of powder): at least 0.5 to 1 cubic cm of the powder should be provided. (The sample will be put into a special quartz capillary (1 or 2 mm in OD, 0.1 mm in wall thickness, and 80 mm in length, which is imported from USA), and its one end will be sealed with a tape.)
- 固態樣品需無毒、無(或極低)揮發性、且無輻射性 (The sample should be free from poison, volatility, and radiation.)  
室溫測量條件下(20°C)，低黏度液體樣品溶液： at least 0.5 to 1 cubic cm of the low-viscosity liquid sample should be provided. (The sample will be put into a special quartz capillary (1 or 2 mm in OD, 0.1 mm in wall thickness, and 80 mm in length, which is imported from USA), and its one end will be sealed with epoxy.) Please be noted that the user should also provide at least 0.5 to 1 cubic cm of pure solvent sample for the X-ray scattering experiment in order to facilitate the background subtraction from the sample solution.
- 樣品若需沈積在基板上，基板的材質及厚度，不可讓 X-ray 對基板的吸收度高於 90% 以上，否則樣品將無法在穿透模式下測量。建議使用 0.1mm 厚度以下的矽晶片或石英玻璃片為佳。(If the sample has to be deposited on a substrate, the material and thickness of the substrate should not cause the degree of absorption of the incident X-ray greater than 90%. Otherwise, the sample will probably not be able to be measured its scattered intensity under the transmission mode of the test. It is suggested to use silicon wafer or quartz glass plate with a thickness less than 0.1 mm.)

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Service Time for High Resolution SAXS System

	Mon	Tue	Wed	Thur	Fri
8:30 am-12:30 noon	3	2	1	2	1
13:30 pm-17:30 pm	3	2	1	2	1

- 校內貴儀使用者預約 (reserved for faculty members or graduate student at NTUST having **MOST** account entitled to use precious instruments of the **MOST**-sponsored labs)
- 校外貴儀使用者預約 (reserved for faculty members or graduate student from other university or research organization having **MOST** account entitled to use precious instruments of the **MOST**-sponsored labs)
- 願意繳交現金之使用者預約 (reserved for applicants willing to pay cash to use the SAXS system)
- 每月第二及第四個週一，儀器停機一天，更換針孔系統及光束停止點。遇假日則順延至次一個上班日。(For the second and fourth Monday of each month, the SAXS system will be down for changing the pinhole and beamstop. Should it be a holiday for that Monday, the down time for the SAXS system will be delayed to the nearest working day.)

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